




МеталлДизайн



2012

Laboratory furniture:
EUROMAX[®], EUROMINI[®], INOVA[®]



МеталлДизайн



For 10 years now specialists of the laboratory furniture factory "METALLDESIGN" study specific features and implement projects of providing laboratory furniture to laboratories of fuel and energy complex, metallurgy, medicine, food, environment and education.

We design furniture, creating new models focused on different needs of enterprises in Russia and countries near and far abroad.

METALDESIGN..... 6

EuroMax..... 9

Side Tables.....	10
Technological Racks.....	14
Titration Rack.....	18
Island Tables.....	20
Technological Racks for Island Tables.....	24
Pedestals.....	30
Cocks and Plumbing Equipment for Side and Island Racks.....	33
Accessories.....	34
Tables with Washers.....	36
Splash Screens.....	40
Table for Sample Preparation with Ventilation.....	42
Mobile Instrument Table.....	44
Weighing Tables.....	46

EuroMini..... 49

Side Tables.....	50
Attached Table.....	52
Shelving Racks.....	54
Titration Rack.....	57
Electrical Equipment for Side Tables.....	58
Accessories.....	60
Pedestals.....	62
Instrument Tables.....	64
Tables with Washers.....	66
Additional Tables.....	70
Weighing Tables.....	74

INOVA.....77

Components for Side and Island Tables.....	78
Accessories.....	84
Order Examples: Lab Side Tables.....	86
Order Examples: Lab Island Tables.....	88
Components for End and Corner Tables.....	90
Order Examples: End and Corner Tables.....	92
Components for Tables with Washers.....	93
Order Examples: Tables with Washers.....	95
Components for End and Corner Tables with Washers.....	96
Order Examples: End Tables with Washers.....	98
Order Examples: Corner Tables with Washers.....	99
Dryers.....	100
Components for Desks.....	101
Order Examples: Desks.....	103
Roller Pedestals.....	104
Shelving Units.....	105
Chairs and Stools.....	106



THOUGHT

Research and development

Since the foundation of the laboratory furniture factory "METALLDESIGN", its specialists study features of using the laboratory furniture and fume hoods in laboratories of fuel and energy complex, metallurgy, medicine, food, environment and education.




CONCEPT

Design and Test

We design furniture, creating new models focused on different needs of enterprises in Russia and countries near and far abroad.

"METALLDESIGN" factory keeps experimental production workshops, where new products are designed, manufactured and tested.

We also test new furniture samples in sectoral research institutes and leading enterprises of Russia.



REALITY

Serial production and modern technology

We use the most advanced software and the most advanced technologies for the production of metal laboratory furniture. We carefully monitor the quality of products at each stage of production.




QUALITY


Due to numerous tests and trials, we can monitor the quality of our products and work on improving it.

The Quality Management System of the factory is certified according to international standard ISO 9001:2008 and accredited in the US, UK, Germany and France.

DESIGN PROCESS



We make the interior design of a future laboratory as two-dimensional drawings and three-dimensional graphics, which allows the customer to consider all the features of ergonomics, functionality, comfort, to optimize the cost of his order and to take aboard the prospects of the enterprise progress.



FURNITURE ASSEMBLY AT THE CUSTOMER'S SITE

Assembly and connection of the furniture and fume hoods to utilities must be performed by qualified personnel only.


We offer services for assembly of laboratory furniture in the customer's site; as well as perform training of certified specialists.




Table Tops



Laminated plate :
Resistant to dilute acids, alkalis and organic solvents for short periods of time.
Withstands short-term temperature exposure up to 170° C.



Laminated plastic LabGrade:
Chemically stable, non-hygroscopic .
Withstands temperature exposure up to 200° C.



Epoxy composite Durcon:
Particularly stable when in contact with most chemicals, acids, alkalis, solvents. Non-hygroscopic.
Does not support combustion.
Withstands temperature exposure up to 300° C.




Porcelain tiles (artificial stone):
Manufactured on the basis of quartz and epoxy resin.
Particularly stable when in contact with most chemicals.
Particularly stable under prolonged thermal contact.
Does not support combustion.
Withstands temperature exposure up to 300° C.



WASHERS AND SINKS



Stainless steel:
Comes complete with LabGrade tabletop



Durcon:
Comes complete with Durcon and LabGrade tabletops.



Durconsink:
Comes complete with Durcon and LabGrade tabletops.



Polypropylene sink:
Comes complete with Laminate, Ceramic, Porcelain, LabGrade tabletops.



When designing EuroMax furniture series, we hoped that every component would comply not only with traditional requirements of Russian laboratories, but also with various international standards. This is why we assert that EuroMax is the laboratory furniture for the most exacting users.

The following requirements were used in the design:
BSEN 13150, BS 3202, DIN 12916,12918,12920

Lab Side Tables





Basic assembly for chemical lab tables in this series is the following:

- Rugged metal frame, powder coated;
- Height of 900 and 750 mm;
- Six different materials of table tops;
- Maximum table load up to 300 kg;
- Adjustable supports and removable front screens for access to services;
- Multisectional tables available;
- Capability of mounting suspended pedestals, technological racks with shelves and power supply, cocks for gas or water with a sink, and a titration rack.

Lab Side Tables

Lower screens hide water, gas, electric service lines. Screens are opened with plastic keys (universal for all side and island tables) and can be easily removed, providing access for servicing



750 mm height

Item Name	Table top	Dimensions, mm	Code
Side table H=750mm L=1202mm	Grey Laminate	1202x750x750	311280
	Ceramic	1202x750x750	311230
	LabGrade	1202x750x750	311290
	Durcon	1202x750x750	311240
	Durcon with collar	1216x764x750	311250
	Ceramic granite	1202x750x750	311260
Side table H=750mm L=1402mm	Grey Laminate	1402x750x750	311480
	Ceramic	1402x750x750	311430
	LabGrade	1402x750x750	311490
	Durcon	1402x750x750	311440
	Durcon with collar	1416x764x750	311450
	Ceramic granite	1402x750x750	311460


663

900 mm height


Item Name	Table top	Dimensions, mm	Code
Side table H=750mm L=1202mm	Grey Laminate	1202x750x900	311180
	Ceramic	1202x750x900	311130
	LabGrade	1202x750x900	311190
	Durcon	1202x750x900	311140
	Durcon with collar	1216x764x900	311150
	Ceramic granite	1202x750x900	311160
Side table H=750mm L=1402mm	Grey Laminate	1402x750x900	311380
	Ceramic	1402x750x900	311330
	LabGrade	1402x750x900	311390
	Durcon	1402x750x900	311340
	Durcon with collar	1416x764x900	311350
	Ceramic granite	1402x750x900	311360

Lab Side Tables

Two-sectionside table



All side and island tables are equipped with height adjustable supports, allowing to set sheerly horizontal working surface.



750 mm height


Item Name	Table top	Dimensions, mm	Code
Two-sectionside table H=750 mm L=2374 mm	Grey Laminate	2374x750x750	311282
	LabGrade	2374x750x750	311292
	Durcon	2374x750x750	311242
	Durcon with collar	2388x764x750	311252
	Ceramic granite	2374x750x750	311262
Two-sectionside table H=750 mm L=2774 mm	Grey Laminate	2774x750x750	311482
	LabGrade	2774x750x750	311492
	Durcon	2774x750x750	311442
	Durcon with collar	2788x764x750	311452
	Ceramic granite	2774x750x750	311462

900 mm height


Item Name	Table top	Dimensions, mm	Code
Two-sectionside table H=900mm L=2374mm	Grey Laminate	2374x750x900	311182
	LabGrade	2374x750x900	311192
	Durcon	2374x750x900	311142
	Durcon with collar	2388x764x900	311152
	Ceramic granite	2374x750x900	311162
Two-sectionside table H=900mm L=2774mm	Grey Laminate	2774x750x900	311382
	LabGrade	2774x750x900	311392
	Durcon	2774x750x900	311342
	Durcon with collar	2788x764x900	311352
	Ceramic granite	2774x750x900	311362

Lab Side Tables

Three-sectionside table



All side and island tables have points for grounding. That allows to operate them safely with any devices.



750 mm height

Item Name	Table top	Dimensions, mm	Code
Three-sectionside table H=750mm L=3546mm	Grey Laminate	3546x750x750	311283
	LabGrade	3546x750x750	311293
	Durcon	3546x750x750	311243
	Durcon with collar	3560x764x750	311253
	Ceramic granite	3546x750x750	311263
Three-sectionside table H=750mm L=4146mm	Grey Laminate	4146x750x750	311483
	LabGrade	4146x750x750	311493
	Durcon	4146x750x750	311443
	Durcon with collar	4160x764x750	311453
	Ceramic granite	4146x750x750	311463

900 mm height


Item Name	Table top	Dimensions, mm	Code
Three-sectionside table H=900mm L=3546mm	Grey Laminate	3546x750x900	311183
	LabGrade	3546x750x900	311193
	Durcon	3546x750x900	311143
	Durcon with collar	3560x764x900	311153
	Ceramic granite	3546x750x900	311163
Three-sectionside table H=900mm L=4146mm	Grey Laminate	4146x750x900	311383
	LabGrade	4146x750x900	311393
	Durcon	4146x750x900	311343
	Durcon with collar	4160x764x900	311353
	Ceramic granite	4146x750x900	311363

Technological racks




- Rugged metal frame, powder coated;
- Shelves for glassware made of tempered glass with aluminum profile;
- Capability of rearrangement of shelves (7 levels for tables with height of 900mm and 9 levels for tables with height of 750mm);
- Cable duct with sockets and automatic shutdown, fluorescent lights;
- Metal duct for installation of cocks for water, gas and air.

Technological racks



Item Name	Dimensions, mm	Code
Technological rack for a side table H=750mm	1200x398x1992	342000
	1400x398x1992	342200
Technological rack for a side table H=900mm	1200x398x1992	342100
	1400x398x1992	342300




Note!


- A technological rack can hold up to three cocks.
One cock for water with a drainage sink and two cocks, gas or air, in various combinations.
- Technological racks of appropriate height only can be mounted on tables with height of 750mm and 900mm.
- Without technological racks only service racks, pedestals and titration racks can be mounted on tables.
- Cocks can be set only with the use of technological racks.

Technological racks

Multisection technological racks



Item Name	Table height	Dimensions, mm	Code
Technological rack for two-sectionside table	H=750mm	2370x398x1992	342002
		2770x398x1992	342202
	H=900mm	2370x398x1992	342102
		2770x398x1992	342302
Technological rack for three-sectionside table	H=750mm	3540x398x1992	342003
		4140x398x1992	342203
	H=900mm	3540x398x1992	342103
		4140x398x1992	342303



Note!

- A technological rack can hold up to three cocks. One cock for water with a drainage sink and two cocks, gas or air, in various combinations.
- Technological racks of appropriate height only can be mounted on tables of 750mm and 900mm height.
- Without technological racks only service racks and pedestals can be mounted on tables.
- Cocks can be set only with the use of technological racks.



Laboratory of Tyumen State University
Operation of "METALLDESIGN" fume hoods since 2004




Knauf Laboratory
Usage of furniture "METALLDESIGN" since 2003

Titration rack





- Metal frame, powder coated;
- Shelf for reagents;
- Screen with frosted glass, illuminated with fluorescent lamps;
- Five stands of stainless steel with ten polypropylene clamps for burettes;
- Duct with two sockets, automatic power and illumination shutdown.

Titration rack



Titration rack complete with a side table is designed to perform various analyses by direct or back titration.



Item Name	Dimensions, mm	Code
Titration rack	1200x420x1992	342900


Island tables





- Rugged metal frame, powder coated, with height of 900 mm and 750 mm;
- Five different materials of table tops;
- Maximum load on one table section up to 500 kg;
- Adjustable supports and removable front screens for access to services;
- Multi sectional tables available;
- Capability of mounting suspended pedestals, technological racks, racks with electric sockets, cocks for gas or water with sink.

Island tables

Island tables



All island tables have special plugs, located in the section supports (both right and left), designed to supply the necessary services.



750 mm height


Item Name	Table top	Dimensions, mm	Code
Island table H=750mm L=1202mm	Grey Laminate	1202x1500x750	316280
	LabGrade	1202x1500x750	316290
	Durcon	1202x1500x750	316240
	Durcon with collar	1216x1514x750	316250
	Ceramic granite	1202x1500x750	316260
Island table H=750mm L=1402mm	Grey Laminate	1402x1500x750	316480
	LabGrade	1402x1500x750	316490
	Durcon	1402x1500x750	316440
	Durcon with collar	1416x1514x750	316450
	Ceramic granite	1402x1500x750	316460

900 mm height



Item Name	Table top	Dimensions, mm	Code
Island table H=900mm L=1202mm	Grey Laminate	1202x1500x900	316180
	LabGrade	1202x1500x900	316190
	Durcon	1202x1500x900	316140
	Durcon with collar	1216x1514x900	316150
	Ceramic granite	1202x1500x900	316160
Island table H=900mm L=1402mm	Grey Laminate	1402x1500x900	316380
	LabGrade	1402x1500x900	316390
	Durcon	1402x1500x900	316340
	Durcon with collar	1416x1514x900	316350
	Ceramic granite	1402x1500x900	316360

Island tables

Two-section island table



Island table tops consist of 2,4 or 6 pieces (depending on the amount of sections), which allows to replace only the damaged piece in case of destruction of the working surface.



750 mm height


Item Name	Table top	Dimensions, mm	Code
Two-section island table H=750mm L=2374mm	Grey Laminate	2374x1500x750	316282
	LabGrade	2374x1500x750	316292
	Durcon	2374x1500x750	316242
	Durcon with collar	2388x1514x750	316252
	Ceramic granite	2374x1500x750	316262
Two-section island table H=750mm L=2774mm	Grey Laminate	2774x1500x750	316482
	LabGrade	2774x1500x750	316492
	Durcon	2774x1500x750	316442
	Durcon with collar	2788x1514x750	316452
	Ceramic granite	2774x1500x750	316462

900 mm height


Item Name	Table top	Dimensions, mm	Code
Two-section island table H=900mm L=2374mm	Grey Laminate	2374x1500x900	316182
	LabGrade	2374x1500x900	316192
	Durcon	2374x1500x900	316142
	Durcon with collar	2388x1514x900	316152
	Ceramic granite	2374x1500x900	316162
Two-section island table H=900mm L=2774mm	Grey Laminate	2774x1500x900	316382
	LabGrade	2774x1500x900	316392
	Durcon	2774x1500x900	316342
	Durcon with collar	2788x1514x900	316352
	Ceramic granite	2774x1500x900	316362

Island tables

Three-section island table



We are not limited by manufacturing three-section tables. The specialists of our factory can design and produce tables with four, five or more sections and technological racks for them.



750 mm height

Item Name	Table top	Dimensions, mm	Code
Three-section island table H=750mm L=3546mm	Grey Laminate	3546x1500x750	316283
	LabGrade	3546x1500x750	316293
	Durcon	3546x1500x750	316243
	Durcon with collar	3560x1514x750	316253
	Ceramic granite	3546x1500x750	316263
Three-section island table H=750mm L=4146mm	Grey Laminate	4146x1500x750	316483
	LabGrade	4146x1500x750	316493
	Durcon	4146x1500x750	316443
	Durcon with collar	4160x1514x750	316453
	Ceramic granite	4146x1500x750	316463

900 mm height


Item Name	Table top	Dimensions, mm	Code
Three-section island table H=900mm L=3546mm	Grey Laminate	3546x1500x900	316183
	LabGrade	3546x1500x900	316193
	Durcon	3546x1500x900	316143
	Durcon с бортиком	3560x1514x900	316153
	Ceramic granite	3546x1500x900	316163
Three-section island table H=900mm L=4146mm	Grey Laminate	4146x1500x900	316383
	LabGrade	4146x1500x900	316393
	Durcon	4146x1500x900	316343
	Durcon with collar	4160x1514x900	316353
	Ceramic granite	4146x1500x900	316363

Technological racks for island tables





- Rugged metal frame, powder coated;
- Shelves for glassware made of tempered glass with aluminum profile;
- Capability of the rearrangement of shelves (7 levels-for table with height of 900mm and 9 levels for table with height of 750mm);
- Cable duct with sockets and automatic shutdown, fluorescent lights;
- Metal duct for setting of cocks for water, gas and air.

Technological racks for island tables



Technological rack for an island table

Item Name	Dimensions, mm	Code
Technological rack for an island table H=750mm	1200x712x1252	342700
	1400x712x1252	342800
Technological rack for an island table H=900mm	1200x712x1095	342500
	1400x712x1095	342600





Note!





A technological rack can hold up to three cocks at each side: one cock for water with a drainage sink and two cocks for gas or air, in various combinations. Only two cocks can be set on a technological rack (one at each side).

Technological racks for island tables

Technological racks for multi-section island tables



Item Name	Table height	Dimensions, mm	Code
Technological rack for two-section island table	H=750mm	2370x712x1252	342702
		2770x712x1252	342802
	H=900mm	2370x712x1095	342502
		2770x712x1095	342602
Technological rack for three-section island table	H=750mm	3540x712x1252	342703
		4140x712x1252	342803
	H=900mm	3540x712x1095	342503
		4140x712x1095	342603



Note!

A technological rack of one section can hold up to three cocks at each side: one cock for water with a drainage sink and two cocks for gas or air, in various combinations. Only two cocks can be set on each section of a technological rack (one at each side).



Laboratory of "Sintez-Kauchuk" factory
Usage of "METALLDESIGN" furniture since 2006



End tables



- Rugged metal frame, powder coated;
- Five different materials of table tops;
- Maximum table top load up to 200 kg;
- Adjustable supports and removable front screens for access to services;
- Capability of usage as part of the island tables, and as standalone product.

End tables

End table





Three-door end table



Item Name	Table top	Dimensions, mm	Code
End table	Grey Laminate	1500x800x900	317580
	LabGrade	1500x800x900	317590
	Durcon	1500x800x900	317540
	Durcon with collar	1514x807x900	317550
	Ceramic granite	1500x800x900	317560
Three-door end table	Grey Laminate	1500x800x900	317780
	LabGrade	1500x800x900	317790
	Durcon	1500x800x900	317740
	Durcon with collar	1514x807x900	317750
	Ceramic granite	1500x800x900	317760

Instruments and accessories can be stored not only in the central compartment of the end table, but also in the side aprons.

Standard end table



Item Name	Table top	Dimensions, mm	Code
Standard End table L=750	Grey Laminate	1500x800x750	317280
	LabGrade	1500x800x750	317290
	Durcon	1500x800x750	317240
	Durcon with collar	1514x807x750	317250
	Ceramic granite	1500x800x750	317260
Standard end table L= 900	Grey Laminate	1500x800x900	317180
	LabGrade	1500x800x900	317190
	Durcon	1500x800x900	317140
	Durcon with collar	1514x807x900	317150
	Ceramic granite	1500x800x900	317160

Note! Only a pedestal with three drawers of EuroMax series can be suspended to a table with height of 750mm.

Pedestals




- Rugged metal frame , powder painted;
- Pedestal doors are filled with sound-absorbing material;
- Capability of mounting pedestals in tables in any combination;
- Individual pedestals for lengthy chemical glassware;
- Plastic inserts in drawers with compartments for small items.

Pedestals



Reservoir pedestal and a pedestal with 4 drawers "LUX" feature fully draw-out guides.

Pedestals with drawers can not only be suspended on EuroMax tables series, but can also be used as a standalone product, due to the rollers and the upper lids with handles.



Reservoir pedestal



Pedestal with a shelf



Pedestal with 3 drawers



Pedestal with 4 drawers




Item Name	Dimensions, mm	Code
Reservoir pedestal, left door	500x520x600	332000
Reservoir pedestal, right door	500x520x600	332100
Pedestal with a shelf, left door	500x520x600	332400
Pedestal with a shelf, right door	500x520x600	332300
Pedestal with 3 drawers (rolling-suspended)	500x520x450(550)	332200
Pedestal with 4 drawers (rolling- suspended)	500x520x600(700)	332900
Pedestal with 4 drawers LUX (rolling- suspended)	500x520x600(700)	332950

Note!

Pedestals with height of 600mm are designed to be suspended on tables with height of 900mm. Two-section pedestals are designed to be suspended on tables of corresponding length.



Pedestals

Two-section pedestal with four drawers




Two-section pedestals are equipped with fully draw-out guides.



Maximum load for each drawer 20kg.



Two-section pedestal with a shelf



All pedestals with EuroMax door series are equipped with reinforced door butts.



Item Name	Dimensions, mm	Code
Тумба двухсекционная с 4-мя ящиками	1200x520x600	332750
Тумба двухсекционная с 4-мя ящиками	1400x520x600	332850
Тумба двухсекционная с полкой, две двери	1200x520x600	332700
Тумба двухсекционная с полкой, две двери	1400x520x600	332800
Тумба двухсекционная без полки, две двери	1200x520x600	332500
Тумба двухсекционная без полки, две двери	1400x520x600	332600

Cocks and plumbing equipment for side and island racks

Item Name			Code
Water cock for side racks Durcon sink Syphon Sewage pipe with length of 500 mm			348240
Water cock for side racks Polypropylene sink Syphon Sewage pipe with length of 500 mm			348200
Water cock for island racks Durcon sink Syphon Sewage pipe of 500 mm length			348241
Water cock for island racks Polypropylene sink Syphon Sewage pipe of 500 mm length			348201
Gas cock			348100
Compressed air cock			348300

Upon customer's request, cocks for other gases can be supplied.

Accessories

Wall dryer

The dryer can be secured on a wall or on technological racks of island tables. There is a connector on the lower part of a frame for connecting the hose to drain water drippings from glassware.

Material of frame and holders- ABS plastic.

Rear wall material - steel, powder painted.

Number of holders - 25 pcs.

Maximum load on holders -1 kg.

Product weight - 5 kg.

Note! When using a dryer avoid contact of frame and holders with organic solvents.

Item Name	Dimensions, mm	Code
Wall dryer	435x155x670	146100



Polypropylene tray for reagents

Item Name	Dimensions, mm	Code
Polypropylene tray for reagents, 10 pcs pack	410x205x45	930000



Eye washer



Item Name	Code
Eye washer	340100

Service stand, all-purpose



Item Name	Dimensions, mm	Code
Service stand	220x200x140	642900



Institute of Hydrocarbon Processing Laboratory
Usage of "METALLDESIGN" furniture since 2005

Technological racks with washers



- Rugged metal frame, powder painted;
- Two different table top materials;
- Different arrangement of a washer on table top: right, left, double;
- Maximum table load up to 200kg;
- Adjustable supports;
- Capability of the splash screen mounting;
- Capability of eye washers mounting;
- Availability of pedestals;
- End tables with table top dimensions 1500x800;
- The service compartment takes the entire volume under the washer.

Tables with washers

End table with washer





Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1500x800x900	317692
LabGrade, Durcon sink	1500x800x900	317694
Durcon, Durcon sink	1500x800x900	317644
Durcon with collar, Durconsink	1514x807x900	317654



The service compartment takes the entire volume under the washer, thereby providing excellent access to services. Storage compartments are located on the right and left of the door, in aprons.

Three-door end table with washer




Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1500x800x900	317892
LabGrade, Durcon sink	1500x800x900	317894
Durcon, Durcon sink	1500x800x900	317844
Durcon with collar, Durcon sink	1514x807x900	317854

Note!
Stainless steel sink is not designed for works with acids.

Tables with washers

Table with washer (rightward)






Table top	Dimensions, mm	Code
LabGrade, HC sink	1202x750x900	315392
	1402x750x900	315792
LabGrade, Durcon sink	1202x750x900	315394
	1402x750x900	315794
Durcon, Durcon sink	1202x750x900	315344
	1402x750x900	315744
Durcon with collar, Durcon sink	1216x764x900	315354
	1416x764x900	315754

Note!

Stainless steel sink is not designed for works with acids.

Table with washer (leftward)






Table top	Dimensions, mm	Code
LabGrade, HC sink	1202x750x900	315292
	1402x750x900	315692
LabGrade, Durcon sink	1202x750x900	315294
	1402x750x900	315694
Durcon, Durcon sink	1202x750x900	315244
	1402x750x900	315644
Durcon with collar, Durcon sink	1216x764x900	315254
	1416x764x900	315654

Note!

The service compartment takes the entire volume under table tops in all tables with washers, shelves are unavailable.

Tables with washers

Table with a double washer






Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1202x750x900	315392
	1402x750x900	315792
LabGrade, Durcon sink	1202x750x900	315394
	1402x750x900	315794
Durcon, Durcon sink	1202x750x900	315344
	1402x750x900	315744
Durcon with collar, Durcon sink	1216x764x900	315354
	1416x764x900	315754

Note!

Stainless steel sink is not designed for works with acids.

Table with a double washer and two mixer taps








Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1202x750x900	315292
	1402x750x900	315692
LabGrade, Durcon sink	1202x750x900	315294
	1402x750x900	315694
Durcon, Durcon sink	1202x750x900	315244
	1402x750x900	315644
Durcon with collar, Durcon sink	1216x764x900	315254
	1416x764x900	315654

Note!

The service compartment takes the entire volume under table tops in all tables with washers, shelves are unavailable.



Splash screens

Screen for tables of 1200 mm length



Item Name	Dimensions, mm	Code
Screen for tables 1200x750x900	1200x60x1970	346100

Screen for tables of 1400 mm length



Item Name	Dimensions, mm	Code
Screen for tables 1400x750x900	1400x60x1970	346300

The screen consists of two vertical supports and the box with stiffeners. The box comprises 8 holes with plugs for mounting the wall dryers. "METALDESIGN" factory produces screens for side tables of 750mm height as well.

Item Name	Dimensions, mm	Code
Screen for tables 1200x750x750	1200x60x1970	346000
Screen for tables 1400x750x750	1400x60x1970	346200




Table for sample preparation with ventilation



- Rugged metal frame, powder coated;
- Two table tops, stainless steel, 1.5 mm thick;
- Rack with a shelf and fluorescent lamp;
- Small-sized exhaust system;
- Opening in a table top for removal of samples and dust residues captured by the cyclone, closed by a lid;
- Cyclone with a removable container to collect residues of dust and samples;
- Maximum table load up to 300 kg;
- Adjustable supports

Table for sample preparation with ventilation






Table is designed for preparation of dry samples.

Movable hood sleeve ensures comfortable work conditions for laboratory assistants.

Powerful fluorescent lights perfectly illuminate the work area.






Item Name	Table top	Dimensions, mm	Code
Table for sample preparation with ventilation	Stainless steel	1500x924x1465	310020

Mobile instrument table



- Rugged metal frame, powder coated;
- Six variants of the working surface;
- Maximum table load up to 200 kg;
- Wheel diameter of 100 mm, front wheels with a stopper;
- The bottom table top is of LabGrade material.

Mobile instrument table



750 mm height

Item Name	Table top	Dimensions, mm	Code
Mobile instrument table H=750mm L=1202mm	Grey laminate/LabGrade	1202x750x750	312280
	Ceramic/LabGrade	1202x750x750	312230
	LabGrade/LabGrade	1202x750x750	312290
	Durcon/LabGrade	1202x750x750	312240
	Durcon with collar/LabGrade	1202x750x750	312250
	Ceramic granite/LabGrade	1202x750x750	312260
Mobile instrument table H=750mm L=1402mm	Grey laminate/LabGrade	1402x750x750	312480
	Ceramic/LabGrade	1402x750x750	312430
	LabGrade/LabGrade	1402x750x750	312490
	Durcon/LabGrade	1402x750x750	312440
	Durcon with collar/LabGrade	1402x750x750	312450
	Ceramic granite/LabGrade	1402x750x750	312460

900 mm height

Item Name	Table top	Dimensions, mm	Code
Mobile instrument table H=900mm L=1202mm	Grey laminate/LabGrade	1202x750x900	312180
	Ceramic/LabGrade	1202x750x900	312130
	LabGrade/LabGrade	1202x750x900	312190
	Durcon/LabGrade	1202x750x900	312140
	Durcon with collar/LabGrade	1202x750x900	312150
	Ceramic granite/LabGrade	1202x750x900	312160
Mobile instrument table H=900mm L=1402mm	Grey laminate/LabGrade	1402x750x900	312380
	Ceramic/LabGrade	1402x750x900	312330
	LabGrade/LabGrade	1402x750x900	312390
	Durcon/LabGrade	1402x750x900	312340
	Durcon with collar/LabGrade	1402x750x900	312350
	Ceramic granite/LabGrade	1402x750x900	312360




Weighing tables



- Very rigid metal frame of thick-walled steel pipe;
- Working plate made of natural granite;
- Four rubber dampers with steel cores for securing the working plate;
- Adjustable supports;
- The possibility of building a weighing complex on the basis of small and large weight tables.




Weighing tables

Small weighing table



Item Name	Table top	Code
Small weighing table H=750 mm	granite	214200
Small weighing table H=900 mm	granite	314200

Large weighing table (without a small one)



Item Name	Table top	Dimensions, mm	Code
Large weighing table H=750 mm	Grey laminate	1202x750x750	314280
		1402x750x750	314480
	LabGrade	1202x750x750	314290
		1402x750x750	314490
Large weighing table H=900 mm	Grey laminate	1202x750x900	314180
		1402x750x900	314380
	LabGrade	1202x750x900	314190
		1402x750x900	314390



LLC "Gipronikel Institute"
Usage of "METALLDESIGN" furniture since 2004



Composition of accessibility, mobility and European quality

The basis of our design of tables comprises the opportunity to form a variety of configurations to successfully complement the laboratory rooms with any area and ceiling height. EuroMini Furniture can supplement or be supplemented with EuroMax furniture series, which is targeted at the most demanding customers. EuroMini furniture is versatile and affordable.

Products meet the standards of GOST 16371-93, GOST 19917-93, GOST 51121-97


Side tables




Basic assembly for chemical lab tables in this series is the following:

- Rugged metal frame, powder coated;
- Six different materials of table tops with dimensions of 1200x600 mm and four variants with dimensions 600x600 mm;
- Maximum table load up to 180 kg;
- Adjustable supports and removable front screens for access to services;
- Capability of mounting pedestals, racks with shelves, ducts with power sockets; cocks for gas or water with sink.

Side tables




Side table is designed for conducting physical or chemical analyses using appropriate instruments (pH tester, spectrophotometer, photocolormeter, etc.)



Item Name	Table top	Dimensions, mm	Code
Side table	Grey laminate	1200x600x900	211180
Side table	Stainless steel	1200x600x900	211120
Side table	Ceramic	1200x600x900	211130
Side table	LabGrade	1200x600x900	211190
Side table	Durcon	1200x600x900	211140
Side table	Ceramic granite	1200x600x900	211160



Attached table



600

600

900



Item Name	Table top	Dimensions, mm	Code
Attached table	Grey laminate	600x600x900	244180
Attached table	LabGrade	600x600x900	244190
Attached table	Durcon	600x600x900	244140
Attached table	Ceramic granite	600x600x900	244160





Laboratory of Analytical Center of water quality control, CJSC "Rosa"
Usage of "METALLDESIGN" furniture since 2005

Racks with shelves




- Rugged metal frame, powder coated;
- Shelves with collars, perfect for storage and accommodation of glassware;
- Removable screens under the shelves.

Racks with shelves




Item Name	Dimensions, mm	Code
Rack with shelves	1180x210x1600	242200
Racks with shelves for attached table	580x210x1600	242000

Racks with sliding glass door shelf




Shelf with sliding doors, 380mm height, provides storage of documents, glassware and non-aggressive chemical reagents.




Item Name	Dimensions, mm	Code
Rack with sliding glass door shelf	1200x360x1600	242300

Titration rack





Item Name	Dimensions, mm	Code
Titration rack	1200x200x1600	242900




- The metal frame, powder coated with a shelf for reagents on top;
- Screen with frosted glass, illuminated with fluorescent lamps;
- Five stands of stainless steel with ten polypropylene clamps for burettes;
- Duct with two sockets, automatic power and illumination shutdown.

Electrical equipment for side tables

Racks with socket blocks




Item Name	Dimensions, mm	Code
Rack with 2 socket block	1200x120x450	242102
Rack with 4 socket block	1200x120x450	242104



Electrical equipment for side tables

Duct with sockets



Item Name	Dimensions, mm	Code
Duct with 2 sockets, 220 V, 1 automatic circuit breaker	1185x90x123	247100
Duct with 4 sockets, 220 V, 1 automatic circuit breaker	1185x90x123	247200

Lamp with a clamp-on

Item Name	Code
Lamp with a clamp-on	247300



Accessories

Side tables accessories



Accessories

Cocks for gas and water

Item Name	Code
Cock for gas with flexible connection	248100
Cock for compressed air with flexible connection	248300
Cock for water with flexible connection, polypropylene sink, siphon, adapter	248200

Note!

When installing a water cock, a single pedestal and a two-socket duct can be suspended on a side table only.

Wall dryer

The dryer can be secured on a wall. The lower part of a frame has a connector to which the hose is connected to drain water drippings from glassware.

Frame and holders material - ABS plastic.

Rear wall material - steel, powder painted.

Number of holders - 25 pcs.

Maximum load on holders - 1 kg.

Product weight - 5 kg.

Item Name	Dimensions, mm	Code
Wall dryer	435x155x670	146100

Service stand, all-purpose

Item Name	Dimensions, mm	Code
Service stand	220x200x140	642900

Polypropylene tray for reagents






Item Name	Dimensions, mm	Code
Polypropylene tray for reagents	410x205x45	930000

Eye washer

Item Name	Code
Eye washer	340100

Shelf with sliding glass doors (for wall mounting)

Item Name	Dimensions, mm	Weight, kg	Code
Shelf with sliding glass doors	1200x330x380	29	241100




Pedestals



- Rugged metal frame, powder painted;
- Pedestal door are equipped with a key lock;
- Pedestal doors are filled with sound-absorbing material;
- Capability of mounting pedestals on tables in any combination;
- Capability of mounting pedestals on attached tables

Pedestals



Item Name	Dimensions, mm	Code
Pedestal(right), with door, 1 shelf	520x520x650	232300
Pedestal(left), with door, 1 shelf	520x520x650	232400
Pedestal with 3 drawers	520x520x650	232200
Pedestal with 5 drawers	520x520x650	232700


Instrument tables



- Rugged metal frame, powder coated;
- Six different materials of table tops with dimensions 1200x600 mm;
- Maximum table load up to 300 kg;
- Adjustable supports and removable front screens for access to services;
- Capability of equipping an extra shelf made of laminated chipboard;
- Full compatibility with all racks and accessories of side tables.


Instrument tables

Instrument table



Item Name	Table top	Dimensions, mm	Code
Instrument table	Grey laminate	1200x600x900	212280
Instrument table	Stainless steel	1200x600x900	212220
Instrument table	Ceramic	1200x600x900	212230
Instrument table	LabGrade	1200x600x900	212290
Instrument table	Durcon	1200x600x900	212240
Instrument table	Ceramic granite	1200x600x900	212260

Instrument table with a shelf



Item Name	Table top	Dimensions, mm	Code
Instrument table with a shelf	Grey laminate	1200x600x900	212380
Instrument table with a shelf	Stainless steel	1200x600x900	212320
Instrument table with a shelf	Ceramic	1200x600x900	212330
Instrument table with a shelf	LabGrade	1200x600x900	212390
Instrument table with a shelf	Durcon	1200x600x900	212340
Instrument table with a shelf	Ceramic granite	1200x600x900	212360

Tables with washers



- Rugged metal frame, powder painted;
- Three different table top materials with dimensions of 1200x600 mm;;
- Different arrangement of a washer on table top: right, left, double;
- Maximum table load up to 100kg;
- Adjustable supports;
- Capability of the splash screen equipping.

Tables with washers






Table top	Dimensions, mm	Code
Stainless steel	500x600x900	215120
Fiberglass	500x600x900	215170
Durcon	500x600x900	215140


Note!

Stainless steel sink is not designed for works with acids. Do not use abrasives to clean the fiberglass sink.

Screen for table 500x600x900

Item Name	Dimensions, mm	Code
Screen for table 500x600x900	500x93x1612	246000

The screen is designed for protection of walls against splashes, and for mounting of a wall dryer.



Tables with washers

Table with a single left washer





Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1200x600x900	215292
Fiberglass	1200x600x900	215270
Durcon	1200x600x900	215244

Table with a single right washer





Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1200x600x900	215392
Fiberglass	1200x600x900	215370
Durcon	1200x600x900	215344

Note!

Stainless steel sink is not designed for works with acids. Do not use abrasives to clean the fiberglass sink.

Tables with washers

Table with a double washer





Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1200x600x900	215492
Fiberglass	1200x600x900	215470
Durcon	1200x600x900	215444

Table with a double washer and two cocks





Table top	Dimensions, mm	Code
LabGrade, stainless steel sink	1200x600x900	215592
Fiberglass	1200x600x900	215570
Durcon	1200x600x900	215544


Note!

Stainless steel sink is not designed for works with acids.
Do not use abrasives to clean the fiberglass sink.

Screen for table 1200x600x900

Item Name	Dimensions, mm	Code
Screen for table 1200x600x900	1180x93x1612	246200

The screen is designed for protection of walls against splashes, and for mounting of the wall dryer.




Additional tables




- Rugged metal frame, powder coated;
- Table tops made of laminated chipboard or LabGrade;
- Maximum table load up to 180 kg;
- Maximum load on a movable table up to 100 kg;
- Adjustable support;
- Desks with pedestals, lockable;
- Computer desks with compartments for system units and with sliding shelves for keyboards.

Additional tables

Desk





Desk with a pedestal



Item Name	Table top	Dimensions, mm	Code
Desk	Grey laminate	1200x600x750	243280
Desk with a left pedestal	Grey laminate	1200x600x750	243380
Desk with a right pedestal	Grey laminate	1200x600x750	243480
Desk	LabGrade	1200x600x750	243290
Desk with a left pedestal	LabGrade	1200x600x750	243390
Desk with a right pedestal	LabGrade	1200x600x750	243480

Additional tables


Computer desk



Item Name	Table top	Dimensions, mm	Code
Left Computer desk	Grey laminate	1200x600x750	243180
Right Computer desk	Grey laminate	1200x600x750	243080
Left Computer desk	LabGrade	1200x600x750	243190
Right Computer desk	LabGrade	1200x600x750	243090


Additional tables

Mobile table




Item Name	Table top	Dimensions, mm	Code
Mobile table	Grey laminate	600x600x800	245180

Mobile instrument table



Lower mobile instrument table



Item Name	Table top	Dimensions, mm	Code
Mobile instrument table	Grey laminate/LabGrade	800x600x900	245880
Mobile instrument table	Durcon/LabGrade	800x600x900	245840
Mobile instrument table	LabGrade/LabGrade	800x600x900	245890
Lower mobile instrument table	Grey laminate/LabGrade	800x600x450	245980
Lower mobile instrument table	Durcon/LabGrade	800x600x450	245940
Lower mobile instrument table	LabGrade/LabGrade	800x600x450	245990




Weighing tables



- Very rigid metal frame of thick-walled steel pipe;
- Working plate made of natural granite of 60mm thickness;
- Four rubber dampers with steel cores for securing the working plate;
- Adjustable supports;
- The possibility of building a weighing complex on the basis of small and large weight tables.



Weighing tables

Small weighing table



Item Name	Table top	Dimensions, mm	Code
Small weighing table	Granite	450x630x750	214200

Large weighing table (without a small one)



Item Name	Table top	Dimensions, mm	Code
Large weighing table	Grey laminate	1200x600x750	214180
Large weighing table	LabGrade	1200x600x750	214190



Kirishi oil refinery plant laboratory
Usage of "METALLDESIGN" furniture since 2004

Weighing tables



Laboratory furniture of the INOVA series has a block-modular construction. INOVA laboratory tables have two basic size in length of 1040 mm and 1560 mm and are assembled from the following Components :

- tabletops;
- pedestals;
- complete sets for side or island tables;
- Complete set for side or island sections;
- racks;
- accessories.

Weighing tables

Table tops




Table tops made of laminated chipboard have durable plastic edge of 2mm thick.

Table top	Dimensions, mm	Code
Laminated chipboard, black	1040x750	600511
Durcon black	1040x750	600510
LabGrade	1040x750	600531
Durcon black with an opening for a sink	1040x750	600510-1
LabGrade with an opening for a sink	1040x750	600531-1





Table top	Dimensions, mm	Code
Laminated chipboard, black	1560x750	600512
Durcon black	1560x750	600506
LabGrade	1560x750	600507
Durcon black with an opening for a sink	1560x750	600506-1
LabGrade with an opening for a sink	1560x750	600507-1


Components of side and island tables

Pedestals

Kam. № 600280



Kam. № 600281




Pedestals are the basis of Inova furniture series, table tops and tables components are attached to them.

All pedestals can be locked with a key.

Useful storage capacity is 160 liters.

Kam. № 600250



Kam. № 600270




Pedestals supports are adjusted from 0 to 60mm, which allows to set it early horizontal working surface.

Item Name	Code
Left pedestal with drawers	600280
Right pedestal with drawers	600281
Pedestal with 3 drawers	600250
Pedestal with 5 drawers	600270


Components of side and island tables

Complete sets for side tables

Kam. № 600001



Kam. № 600002




Complete sets are necessary for connecting pedestals and frame tension clamps. Complete set also includes decorative screens to close the technological sections and fittings for tables assembly.


Item Name	Code
Complete set for Inova side tables 1040	600001
Complete set for Inova side tables 1560	600002

Complete sets for additional side sections

Kam. № 600003



Kam. № 600004





Additional side sections are designed for maximizing space for legs at seated activity.

Item Name	Code
Complete set for additional Inova side sections 1040	600003
Complete set for additional Inova side sections 1560	600004

Components of side and island tables



Complete sets for island tables



In addition to the listed above complete sets also include grounding straps, to which a terminal box and a power cord are attached.

Item Name	Code
Complete set for Inova island tables 1040	600005
Complete set for Inova island tables 1560	600006

Complete sets for additional island sections





Decorative panels of island tables are equipped with special plugs of large diameter which can be removed in case of connecting services to island tables.

Item Name	Code
Completes et for additional Inova island sections 1040	600007
Completes et for additional Inova island sections 1560	600008

Components of side and island tables

Side racks with shelves




All racks in the Inova series are equipped with fluorescent lights with screens.

Power is supplied from a terminal box, which is part of the table elements.

Item Name	Dimensions, m	Code
Side rack with a shelf Inova1040	1040x330x530	641100
Side rack with a shelf Inova1560	1560x330x530	641400

Racks with sliding glass door shelves



Racks mounting and fixing is performed with the help of a table top, which already has holes with plugs.

Item Name	Dimensions, m	Code
Rack with sliding glass door shelves Inova1040	1040x350x1100	641200
Rack with sliding glass door shelves Inova1560	1560x330x1100	641500

Components of side and island tables

Rack with metal doors shelf



Item Name	Dimensions, mm	Code
Rack with metal doors shelf Inova	1040x350x530	641300

Service stand, all-purpose



Stand can be mounted anywhere on a table top.
Holes for securing and power supply are made on-the-spot.


Item Name	Dimensions, mm	Code
Service stand	220x200x140	642900

5 power sockets of 220 V, 2 RJ-45 connectors, 16A circuit breaker


Accessoires

Cock for water with a drainage sink and plumbing equipment

Kam. N° 648240



Kam. N° 648200



Item Name	Code
Cock for water with a Durcon sink for Durcon table top	648240
Cock for water with a polypropylene sink for LabGrade table top	648200

Gas cocks



You can order an unlimited number of gas cocks, and place them any where on the work surface.

Note!

Vacuum cock is supplied without a connection!

Item Name	Code
Vacuum cock Gas	648000
Cock with a flexible connection	648100
Compressed air cock with a flexible connection	648300



Henkel-Era Laboratory of Kirishi
Usage of "MetallDesign" furniture since 2004

Order examples: Lab side tables

Side table 1040 (a pedestal with 5 drawers and a pedestal with a door)



The type of pedestals and tabletops can be changed at will in the side and island tables.

Product kit	Code	Pcs
Laminated chipboard, black	600511	1
Pedestal with 5 drawers	600270	1
Pedestal with shelves, right	600281	1
Completes et for Inova side tables 1040	600001	1

Side table 1560 (pedestal with 3 drawers and a pedestal with a door)




Product kit	Code	Pcs
Laminated chipboard, black	600512	1
Pedestal with 3 drawers	600250	1
Pedestal with shelves, right	600280	1
Complete set for Inova side tables 1040	600002	1

Order examples: Lab side tables


Side table 1040 (pedestal with 5 drawers, pedestal with door, rack with a glass shelf)

Product kit	Code
Laminated chipboard, black	600511
Pedestal with 5 drawers	600270
Pedestal with shelves, right	600281
Side rack with a shelf Inova1040	641100
Completes et for Inova side tables 1040	600001



Side table 1560 (pedestal with 3 drawers, pedestal with door, rack with sliding doors)

Product kit	Code
Laminated chipboard, black	600512
Pedestal with 3 drawers	600250
Pedestal with shelves, left	600280
Rack with sliding glass doors shelf Inova1560	641500
Completes et for Inova side tables 1560	600002



Order examples: Lab island tables

Island table 1040 (pedestals with right and left doors)



Product kit	Code	Pcs
Laminated chipboard, black	600511	2
Pedestal with shelves, left	600280	2
Pedestal with shelves, right	600281	2
Completes et for island tables Inova 1040	600005	1

Island table 1560 (pedestals with 5 drawers and pedestals with doors)





Table tops of island tables are joined together out of two side table tops. Table tops can be combined at will, but be careful: when joining LabGrade and Durcon table tops with laminate ones, the ledge of 10mm thick between different table tops is unavoidable.

Product kit	Code	Pcs
Laminated chipboard, black	600512	2
Pedestal with 5 drawers	600270	2
Pedestal with shelves, left	600280	2
Completes et for island tables Inova1560	600006	1


Order examples: Lab island tables

Island table 1040 (pedestals with right and left doors, rack with a glass shelf, rack with metal doors shelf and lights)



Product kit	Code	Pcs
Laminated chip board, black	600511	2
Pedestal with shelves, left	600280	2
Pedestal with shelves, right	600281	2
Side rack with a shelf Inova1040	641100	1
Rack with metal doors shelf Inova	641300	1
Complete set for island tables Inova1040	600005	1

Island table 1560 (pedestals with right and left doors, rack with a glass shelf, rack with sliding glass doors shelf)




Product kit	Code	Pcs
Laminated chip board, black	600512	2
Pedestal with shelves, left	600281	2
Pedestal with shelves, right	600280	2
Side rack with a shelf Inova1560	641400	1
Rack with sliding glass doors shelf Inova1560	641500	1
Complete set for island tables Inova1560	600006	1

Components of end and corner tables

Table tops


End tables

Table top	Dimensions, mm	Code
Laminated chip board, black	1500x750	600545
Durcon black	1500x750	600544
LabGrade	1500x750	600546




Exterior corner tables

Table top	Dimensions, mm	Code
Laminated chip board, black	750x750	600542
Durcon black	750x750	600541
LabGrade	750x750	600543



Interior corner tables

Table top	Dimensions, mm	Code
Laminated chip board, black	970x750	600535
Durcon black	970x750	600534
LabGrade	970x750	600536



Pedestals




Table top	Code
Corner pedestal with a left door	600295
Corner pedestal with a right door	600296


Components of end and corner tables

Complete sets for end and corner tables


Kam. № 600009



Kam. № 600010



Kam. № 600011




Item Name	Code
Complete set for interior corner table	600009
Complete set for exterior corner tables	600010
Complete set for end tables	600011

Order examples: end and corner tables

End table

Product kit	Code	Pcs
Laminated chip board, black	600545	1
Corner pedestal with a left door	600295	1
Corner pedestal with a left door	600296	1
Complete set for end tables	600011	1

When combining island sections of Inova furniture with end tables you get a fully-featured table.



Exterior corner table

Product kit	Code	Pcs
Laminated chip board, black	600542	1
Corner pedestal with a right door	600296	1
Complete set for exterior corner tables	600010	1


Exterior corner tables allow to create a continuous side table, even when there are exterior angles in the laboratory.



Interior corner table

Product kit	Code	Pcs
Laminated chip board, black	600535	1
Corner pedestal with a right door	600296	1
Complete set for interior corner tables	600009	1

Internal corner tables are necessary to fill the corners of laboratories and to create continuous tables.




Components of tables with washers

Table tops

Side table with washer, small

Table top	Dimensions, mm	Code
Durcon black with one washer	520x750	600597


Kam. № 600597




Side table with washer, small

Table top	Dimensions, mm	Code
Durcon black	1040x750	600598

Kam. № 600598



Pedestals



Washer pedestals differ from conventional ones in reinforced frame and the lack of locks.


Item Name	Code
Washer pedestal with a left door	600310
Washer pedestal with a right door	600311
Double washer pedestal	600320




Components of tables with washers

Complete sets for tables with washers


Kam. № 600101



Kam. № 600118



Kam. № 600103



Components of tables with sinks, in addition to screens, fasteners and ties include drain sinks with sanitary equipment and mixer taps with flexible connections.

Drain sinks of Inova series are suspended to the bottom of a tabletop by means of special brackets.

Item Name	Code
Complete sets for small side tables with washers	600101
Complete sets for large side tables with washers	600118
Complete sets for island tables with washers	600103

Order examples: tables with washers

Tables with washers can be used as a standalone product and as part of side or island tables.

Small side table with a washer 520x750x900

Product kit	Code	Pcs
Durcon black 520x750	600597	1
TWasher pedestal with a left door	600310	1
Complete sets for small side tables with washers	600101	1




Large side table with a washer 1040x750x900

Product kit	Code	Pcs
Durcon black 1042x750	600598	1
Double washer pedestal	600320	1
Complete set for large side tables with washers	600118	1



Island table with washers 520x1500x900

Product kit	Code	Pcs
Durcon black 520x750	600597	2
Washer pedestal with a left door	600310	1
Washer pedestal with a right door	600311	1
Complete set for island tables with washers	600103	1




Components of end and corner tables with washers

Table tops


End tables

Working surface	Dimensions, mm	Code
Durcon black with one washer	1500x750	600595
Durcon black with two washers	1500x750	600594

Kam. № 600595




Kam. № 600594



External corner tables with washers

Working surface	Dimensions, mm	Code
Durcon black	750x750	600596


Kam. № 600596



Internal corner tables with washers

Working surface	Dimensions, mm	Code
Durcon black	970x970	600599

Kam. № 600599



Pedestals




Item Name	Code
Corner washer pedestal with a left door	600300
Corner washer pedestal with a right door	600301
Double washer pedestal	600320


Components of end and corner tables with washers

Complete sets for end and corner tables with washers


Kam. № 600121




Kam. № 600124



Kam. № 600122




Kam. № 600107



Item Name	Code
Complete set for end tables with one washer	600121
Complete set for end tables with two washers	600124
Complete set for external corner tables with washers	600122
Complete set for internal corner tables with washers	600107


Order examples: end tables with washers

End table with one washer



Product kit	Code	Pcs
Durcon black with one washer	600595	1
Double washer pedestal	600320	1
Complete set for end tables with one washer	600121	1

End table with two washers



Product kit	Code	Pcs
Durcon black with two washers	600594	1
Corner washer pedestal with a left door	600300	1
Corner washer pedestal with a right door	600301	1
Complete set for end tables with two washers	600124	1


Order examples: corner tables with washers

External corner table with a washer



Product kit	Code	Pcs
Durcon black	600596	1
Corner washer pedestal with a right door	600301	1
Complete set for external corner tables with washers	600122	1


Internal corner table with a washer



Product kit	Code	Pcs
Durcon black	600599	1
Corner washer pedestal with a left door	600300	1
Complete set for internal corner tables with washers	600107	1

Dryers

Table dryer



Item Name	Dimensions, m	Code
Table dryer	520x245x785	646100

Designed for mounting on table tops w/codes
600595, 600597, 600598, 600599

Otherwise wall dryers can be used

Wall dryer

Wall dryer
Frame and holders material - ABS plastic.
Rear wall material - steel, powder painted.
Number of holders - 25 pcs.
Maximum load on holders -1 kg.
Product weight - 5 kg.

Note!

**When using a dryer avoid contact of frame and holders
with organic solvents.**



Item Name	Dimensions, m	Code
Wall dryer	435x155x670	146100

Components for desks

INOVA desks have one basic length dimension equal to 1560 mm, and two height dimensions: 750 and 900 mm.

Table tops

Working surface	Dimensions, mm	Code
Laminated chipboard, black	1560x750	600512


For 750 mm desks

Pedestals


For 900 mm desks




Kam. № 600245




Kam. № 600246



Kam. № 600280



Kam. № 600281




Item name	Table height	Code
Pedestal with shelves, right	750	600246
Pedestal with shelves, left	750	600245
Pedestal with drawers	750	600240
Pedestal with shelves, left	900	600280
Pedestal with shelves, right	900	600281
Pedestal with 3 drawers	900	600250
Pedestal with 5 drawers	900	600270


Components for desks

Complete set for a desk


Kam. № 600012




Kam. № 600015



Kam. № 600013




Kam. № 600017



Item name	Table height	Code
Complete set for a desk (pedestal leftwards)	750	600012
Complete set for a desk (pedestal rightwards)	750	600013
Complete set for a desk (pedestal leftwards)	900	600015
Complete set for a desk (pedestal rightwards)	900	600017

Order examples: desks

Desk with height of 750 mm (pedestal with shelves, left)




Any table tops with length of 1560mm as well as racks and accessories can be mounted on desks.

Note!
A desk with height of 750 mm can not be joined to a desk with height of 900 mm

Product kit	Code
Laminated chipboard, black	600512
Pedestal with shelves, left	600245
Complete set for a desk (pedestal to the right)	600012

Desk with height of 750 mm (pedestal with drawers)



Note!
Roller pedestals are not designed for desks with height of 750 mm

Product kit	Code
Laminated chipboard, black	600512
Pedestal with drawers	600240
Complete set for a desk (pedestal to the left)	600013

Roller pedestals

Item name	Dimensions, mm	Code
Roller pedestal with a shelf, right	520x500x770	632300



Item name	Dimensions, mm	Code
Roller pedestal with a shelf, left	520x500x770	632400



Item name	Dimensions, mm	Code
Roller pedestal with 3 drawers	520x500x770	632200




Item name	Dimensions, mm	Code
Roller pedestal with 4 drawers	520x500x770	632900





Shelving units

Flower shelving unit



Shelving unit



Item name	Shelves	Dimensions	Code
Flower shelving unit	Tempered glass	1000x540x2000	734000
Shelving unit	Metal	1000x540x2000	634000

Chairs and stools


Laboratory chairs and stools



Item Type	Seat cover	Height (min-max)	Code
Chair	Polyutherane	540-670	147700
Chair with rollers	Vinyl	540-670	147100
Chair with supports	Vinyl	540-670	147200
Stool with rollers	Vinyl	560-690	147400
Stool with supports	Vinyl	560-690	147500
Stool with rollers	Polyutherane	560-690	147800
Stool with supports	Polyutherane	560-690	147900

Office chairs (for work at desks with height of 720-750 mm)

Item Type	Seat cover	Height (min-max)	Code
Chair with rollers	Vinyl	410-540	147300
Stool with rollers	Vinyl	410-540	147600



195197, St. Petersburg, Polustrovskiy, 50.
Tel.: 8 (812) 703-01-65
e-mail: info@lenlab.ru
www.labmax.ru



AB24

ПР 019